




Lohmann & Rauscher

SwellSpots™

Oedema and fibrosis pads.



Ideal for difficult areas of oedema, fibrosis, skin folds, misshapened limbs as well as padding and protection.



POWERED BY
coolcore®

Overview

Why choose SwellSpots™?

POWERED BY


SwellSpots, oedema and fibrosis pads are designed to address particularly troublesome areas by:

- softening tissue fibrosis/hardened skin
- filling concavities
- padding bony prominences
- diffusing pressure points
- protecting vulnerable areas

SwellSpots are available in a wide range of versatile shapes and sizes that tuck into compression garments or under bandages. Powered by Coolcore, these SwellSpots rapidly wick away sweat and moisture from the body, keeping it cool and comfortable.

What is Coolcore?

Coolcore is an innovative, chemical-free fabric that rapidly wicks away sweat and body moisture to regulate body temperature. Its chemical-free technology is built into the fibres of the fabric to withstand its cooling properties through wash and wear during the one year life of the garment.




88%
faster wicking
sweat and moisture
away from the body*



55%
faster
drying time
for easier
garment care*



30%
more
cooling power
to better regulate
body temperature*




foam technology helps soften
tissue fibrosis and hardened skin



creates multipoint compression
profiles which help encourage
lymphatic reabsorption




versatile shapes and sizes
available off-the-shelf



the special Coolcore fibre
helps wick away moisture



saves time from cutting
your own chip bags




chevron channelling based on
manual lymphatic drainage (MLD)
principles of directing lymphatic
fluid to collateral pathways

* Tested by Hohenstein Institute, an accredited laboratory and research institute focused on testing and certifying textiles.

SwellSpots™

Head-Neck Pad, Versi-Face Pad, Curved Pad & Crescent Pad



Head-Neck Pad (previously Reynolds)



Versatile pad designed for use on the face and neck by placing the Head-Neck Pad in various locations along the strap.

Pad width (cm)	Pad height (cm)	REF	Shipping unit
16.5	7.5	104 544	1/5

Versi-Face Pad

Versatile pad designed for use around the eyes and under the chin.

Width (cm)	Height (cm)	REF	Shipping unit
29	16	104 542	1/5



Curved Pad


Versatile pad designed for use around the ankle, patella, elbow, eyes and behind the ear.

Width (cm)	Height (cm)	REF	Shipping unit
5	13	104 527	1/5

Crescent Pad

Wedge-shaped pad designed to support lobules and skin folds.

Length (cm)	Height (cm)	Diameter	REF	Shipping unit
51	6.5	9	104 547	1/5
61	6.5	11	104 546	1/5
71	6.5	12	104 545	1/5



SwellSpots™

Carpal Hand Pad, Dorsal Pad, Rectangular-Tapered Pad & Knee-Elbow Shaper Pad

Carpal Hand Pad (previously Jo-Ju)

Pad designed to use at the thenar web space.

Width (cm)	Height (cm)	REF	Shipping unit
4.5	7	104 537	1/5



Dorsal Pad

Pad designed to use on the dorsum of either the hand or foot.


Width (cm)	Height (cm)	REF	Shipping unit
10	14	104 529	1/5
10	19	104 528	1/5



Rectangular-Tapered Pad (previously Katina Small & Large)

Tapered pad designed for use at the ulnar aspect of the forearm, the IT band of the lateral thigh, and the anterior surface of the calf.


Width (cm)	Height (cm)	REF	Shipping unit
11	36	104 534	1/5
18	42	104 533	1/5



Knee-Elbow Shaper Pad (previously LT Knee Shaper)

Versatile pad designed for use on the knee and elbow.

Girth range (cm)	REF	Shipping unit
35-65	104 557	1/5



SwellSpots™

Oval Pad, Ring Pad, Kidney Pad & Tic-Tac Pad



Oval Pad



Versatile pad designed for use inside the hand splints, along the ulnar aspect of the forearm, and the groin.

Width (cm)	Height (cm)	REF	Shipping unit
8	19	104 532	1/5

Ring Pad

Versatile pad designed to fill crevices or lift lobules at the ankle, wrist and abdominal pannuses. The Ring Pad is available with or without a VELCRO® brand hook closure.

Width (cm)	Height (cm)	Diameter	REF	Shipping unit
60	6.5	9	104 558*	1/5
60	6.5	9	104 559	1/5



Kidney Pad (previously Kimbe)


Versatile pad designed for use at the scapula, groin and neck.

Width (cm)	Height (cm)	REF	Shipping unit
14	24	104 538	1/5

Tic-Tac Pad

Versatile pad designed to address larger areas of oedema and fibrosis, such as the back, abdomen and thigh.

Size (cm)	REF	Shipping unit
20 x 20	104 563	1/5
30 x 30	104 562	1/5
40 x 40	104 561	1/5



SwellSpots™

Labia Pad, Female Genital Pad, Scrotum Pad & Circular Pad



Labia Pad

Tapered pad designed for the female genital area that can be worn inside or over underwear.

Width (cm)	Height (cm)	REF	Shipping unit
33	53	104 551	1/5

Female Genital Pad (previously Sweet Spot)

Pad designed for use in the female genital area that can be placed inside or over underwear. The black side of the pad is made of non-skid material to help reduce shifting.


Width (cm)	Height (cm)	REF	Shipping unit
7	15	104 560	1/5



Scrotum Pad

Tapered pad designed for the male genital area that can be worn inside or over underwear.

Width (cm)	Height (cm)	REF	Shipping unit
34	53	104 552	1/5



Circular Pad (previously Sobol Small & Large)

Versatile pad designed for use on the chest wall, hip, anterior knee, posterior knee and abdomen.

Diameter (cm)	REF	Shipping unit
18	104 536	1/5
22	104 535	1/5



SwellSpots™

Full Bra Pad, Lateral Bra Pad, Breast Pad & Palmer Pad



Full Bra Pad


Pad designed for use from the lateral chest wall toward the anterior midline of the chest, and can be worn inside a bra.

Cup size*	REF	Shipping unit
A-B cup	104 556	1/5
B-C cup	104 555	1/5
C-D cup	104 554	1/5

Lateral Bra Pad

Pad designed for use on the lateral chest wall, and can be worn inside a bra.

Width (cm)	Height (cm)	REF	Shipping unit
10	16	104 553	1/5



Breast Pad

Pad designed specifically for use in the breast area, and can be worn inside a bra. Its seamless design allows the Breast Pad to be used on either side of the chest.

Cup size*	REF	Shipping unit
A-B cup	104 541	1/5
C-D cup	104 540	1/5
D-DD cup	104 539	1/5

Palmar Pad


Versatile pad designed for use on palmar surface of the hand, dorsum, concave shape of palm and inside a bra.

Width (cm)	Height (cm)	REF	Shipping unit
7	8	104 531	1/5
8	9	104 530	1/5



SwellSpots™


Toe Pad & Thorax-Thigh Pad



Toe Pad (previously We'a Toe Small, Medium & Large)

Unique digit spacers pad designed for comfort and separation of individual toes.

MTP circumference (cm)	REF	Shipping unit
<22	104 550	1/5
<25	104 549	1/5
<32	104 548	1/5



Thorax-Thigh Pad

Versatile pad designed for use at the chest, groin, medial thigh and shoulder.

Width (cm)	Height (cm)	REF	Shipping unit
52	22	104 543	1/5

Care Instructions for SwellSpots

Washing

- machine-wash with mild detergent on permanent press
- do not use fabric softener, bleach or other additives
- do not dry-clean

Drying

- tumble-dry on permanent press
- do not iron

Materials

- blend of: polyurethane, nylon, polyester, elastane
- not made with natural rubber latex



Indications

- lymphoedema (when used with bandages, e.g. Rosidal K/other compression systems, e.g. ReadyWrap)
- oedema (when used with bandages, e.g. Rosidal K/other compression systems, e.g. ReadyWrap)

Precaution

- allergies or sensitivities to any of the product materials
- care should be taken in patients using anticoagulant medication to assess tolerance

For all queries relating to SwellSpots as well as ReadyWrap, Garesia, TributeNight or Exo products please contact your local distributor:

VELCRO® is a registered trademark of Velcro Industries B.V. Used with permission.

Lohmann & Rauscher Pty Ltd · 321 Ferntree Gully Road · Mount Waverley, Victoria 3149, Australia · Customer Service: 1300 572 869
Lohmann & Rauscher Limited · 21 Queen Street · Auckland 1010 · New Zealand · Customer Service: 0800 572 869