



Lohmann & Rauscher

## TributeWrap™ **NEW**

Ease and comfort within your reach.

Self-adjustable compression foam garment to be used at rest or during the night for the management of lymphoedema or oedema.



# TributeWrap™

## Promoting skin health with proven foam technology.

TributeWrap is an adjustable compression garment with foam technology that nestles into the tissue to promote skin health<sup>1</sup>. It supports easy self-management of lymphoedema and oedema during low activity by comfortably adjusting to your patient's unique shape and lifestyle.

TributeWrap helps maintain gains made during therapy, while allowing the flexibility to support further reduction. Simply slide it on, adjust the straps, and cover with Sleep Sleeve for comfortable evening compression.



Off-the-Shelf



Evening & Nighttime



Low Activity

[Explore the other great benefits of TributeWrap >](#)



**Easy-to-use  
finger grip on  
every strap.**

Simple adjustability  
for the right fit  
and contour.




**Proven foam  
technology**  
Breaks up  
and softens  
hardened or  
fibrotic tissue  
to support  
skin health.



**Sleep Sleeve**  
Complimentary  
cover included  
with all styles to  
ensure the straps  
don't get tangled  
in the sheets.

POWERED BY  **coolcore**

Innovative interior fabric that rapidly wicks away sweat and moisture to regulate body temperature with chemical-free Cooling Power.








TributeWrap™ Wrist to Axilla



### Chevron channeling

Designed based on MLD principles of directing fluid to collateral pathways.



### Adjustable strapping system

Easily conforms for the right, comfortable fit, and has flexibility to support further reduction.



### Supports mobility

Angled straps at flexion points for a comfortable and functional design.



### 6-month warranty

TributeWrap is guaranteed to be free of material and manufacturing defects for the six-month life of the garment.



“In the morning, I took it off and couldn’t believe the difference! I wondered if this was a fluke, but I’ve noticed the same results every morning!

This is outstanding! It is easier to put on my daytime flat knit garments, and my legs have less fatigue throughout the day.

I am convinced TributeWrap does make a significant difference!”

John, Lymphoedema patient

# TributeWrap™ NEW

## Self-adjustable compression foam garment.

### Indications

- mild to severe lymphoedema, stage 1–3
- chronic oedema (non-blood related)

### Precautions

- impaired sensation (peripheral neuropathy)
- cognitive impairment
- any open wounds should be dressed prior to application



### Contraindications

- moderate to severe peripheral arterial disease
- arterial insufficiency or ABPI < 0.8
- untreated congestive heart failure
- untreated infection (e.g. cellulitis)
- allergies or sensitivities to materials
- acute DVT
- severe cognitive impairment

### TributeWrap Glove

Self-adjustable compression foam garment for hand and fingers



Size	Circumference (cm)			REF
	Palm (cB)	Wrist (cC)		
S	16 – 20	16.5 – 21	Left	146 783
			Right	146 786
M	18 – 22	17.5 – 22	Left	146 784
			Right	146 787
L	20 – 24.5	18.5 – 23	Left	146 785
			Right	146 788

 Not made with natural rubber latex 

### TributeWrap Arm

Self-adjustable compression foam garment for wrist to axilla


Size	Circumference (cm)			Length (cm)		REF
	Wrist (cC)	Elbow (cE)	Axilla (cG)			
S	16 – 21	22 – 28	25 – 31	43 – 47	Left	146 777
					Right	146 780
M	17 – 22	26 – 32	30 – 37		Left	146 778
					Right	146 781
L	18 – 23	30 – 36	35 – 44		Left	146 779
					Right	146 782

 Not made with natural rubber latex 

### TributeWrap Below Knee

Self-adjustable compression foam garment for foot and calf

Size	Circumference (cm)		Length (cm)			REF
	Ankle (cB)	Calf (cC)	Foot (IA)	Calf (ID)		
S	24 – 30	30 – 36	19 – 22	40 – 44	Left	146 789
					Right	146 792
M	27 – 33	33 – 41.5			Left	146 790
					Right	146 793
L	30 – 36	36 – 45			Left	146 791
					Right	146 794

 Not made with natural rubber latex 