

VENOSAN®

LEGLINE 30

Support Stockings 280 den

ELEGANT LUXURY



VENOSAN®

LEGLINE 30


30 mm Hg
280 den

Support Stockings 280 den

Stützstrümpfe | Support stockings | Bas de support | Calze riposanti



AD
Wadenstrümpfe
Below knee
Bas au mollet
Gambaletti



AGH
Schenkelstrümpfe
Thigh length
Bas à cuisse
Calza alla coscia



AT*
Strumpfhose
Pantyhose/Tights
Collants
Collant


* Komfort-Hosenteil
Comfort panty
Culotte confort
Corpino comfort

Stufenlos abnehmende Kompression | Graduated compression |
Compression dégressive | Compressione graduata

Umfang | Circumference | Circonférence | Circonferenza (cm)

AD	AGH, AT	
	60% 18 mm Hg	cF
80% 24 mm Hg	80% 24 mm Hg	cD
		cC
100% 30 mm Hg	100% 30 mm Hg	cB

XS	S	M	L	XL
42-44	45-47	48-50	51-53	54-60
29-31	32-34	35-37	38-40	41-43
31-33	34-36	37-39	40-42	43-45
20-21	22-23	24-25	26-27	28-29



Effective Compression

Verhindert geschwollene Beine
Prevents swollen legs
Évite des jambes enflées
Impedisce le gambe gonfie

Everyday Elegance

Transparente Optik
Transparent look
Look transparent
Look trasparente

Optimal Comfort

Für höchsten Tragekomfort, rutscht nicht
For maximum comfort, does not slip
Pour un confort optimal, ne glisse pas
Non stringe, non scivola

Spitzenhaftband | Lace top band | Adhésive ornée | Bordo autoreggente pizzo




Aussen
Outside
Face externe
Lato esterno



Innen
Inside
Côte intérieur
Lato interno

Grösse | Height | Grandeur | Altezza



Artikel | Article | Article | Articolo

Type	*Colours	XS	S	M	L	XL
AD	Nude	VL 403 NU	VL 404 NU	VL 405 NU	VL 406 NU	-
	Sahara	VL 403 SA	VL 404 SA	VL 405 SA	VL 406 SA	-
	Black	VL 403 BL	VL 404 BL	VL 405 BL	VL 406 BL	-
AGH	Nude	VL 423 NU	VL 424 NU	VL 425 NU	VL 426 NU	-
	Sahara	VL 423 SA	VL 424 SA	VL 425 SA	VL 426 SA	-
	Black	VL 423 BL	VL 424 BL	VL 425 BL	VL 426 BL	-
AT	Nude	VL 433 NU	VL 434 NU	VL 435 NU	VL 436 NU	VL 437 NU
	Sahara	VL 433 SA	VL 434 SA	VL 435 SA	VL 436 SA	VL 437 SA
	Black	VL 433 BL	VL 434 BL	VL 435 BL	VL 436 BL	VL 437 BL

*Colours



Nude



Sahara



Black

Materialzusammensetzung | Composition | Composition | Composizione

AD|AGH: 70% Polyamide | 30% LYCRA® Elastane

AT: 78% Polyamide | 22% LYCRA® Elastane





VENOSAN®
Compression Stockings

www.venosan.com

VENOSAN® is a registered
trademark of SWISSLASTIC AG
ST. GALLEN, Switzerland

an  Company



SWISSLASTIC AG ST. GALLEN
Rorschacher Strasse 304
CH-9016 St. Gallen
Switzerland